



***Fluids and Combustion Facility  
Preliminary Design Review***

---



# **FCF System**

Chris Pestak  
February 12, 2001




# ***Fluids and Combustion Facility***

## ***Preliminary Design Review***



### **MRDOC Organization**






# Fluids and Combustion Facility

## Preliminary Design Review



### FCF Project Organization under MRDOC





# ***Fluids and Combustion Facility***

## ***Preliminary Design Review***



---

### **FCF Project Organization - FDC Team Contact List**

Christopher J. Pestak	FCF Director	216/925-1050
Andrew M. Peddie	FCF Deputy Director	
	FIR & SAR PDT Lead	216/925-1044
Jon C. Wetherholt	Systems Engineering Lead	216/925-1088
Dr. Nora G. Bozzolo	Common Systems PDT Lead	216/925-1038
Martin A. O'Toole	CIR PDT Lead	216/925-1200
Joel J. Knapp	Logistics & Operations Lead	
	Integration & Utilization Lead	216/925-1060
Stephen M. Lux	Software Lead	216/925-1055
William F. Quinn	Product Assurance Lead	216/925-1293
Dr. Mark W. Pestak	Science Liaison	216/925-1136



***Fluids and Combustion Facility  
Preliminary Design Review***

---



# **FCF Overview**



# ***Fluids and Combustion Facility***

## ***Preliminary Design Review***



### **FCF Acronym List**

AMA	Atmospheric Measurement Assembly	ISPR	International Standard Payload Rack
ARIS	Active Rack Isolation Subsystem	ISS	International Space Station
ATCS	Air Thermal Control System	ITCS	Internal Thermal Control System
ATCU	Air Thermal Control Unit	LED	Light Emitting Diode
CIR	Combustion Integrated Rack	LLL	Low Light Level
DCM	Diagnostics Control Module	LMM	Light Microscopy Module
ECS	Environmental Control Subsystem	MDCA	Multi-user Droplet Combustion Assembly
EDU	Experiment Development Unit	MFC	Mass Flow Controller
EM	Engineering Model	OB	Optics Bench
EMC	Electro- Magnetic Compatibility	OM	Optics Module
EP	Experiment Package	OPI	Optics Plate Interface
EPCU	Electrical Power Control Unit	ORU	Orbital Replacement Unit
EPS	Electrical Power Subsystem	PCS	Portable Computer System
ESP	Electronic Support Package	PDC	Payload Development Center
ESSA	EPCU Shutoff Switch Assembly	PFE	Portable Fire Extinguisher
EVP	Exhaust Vent Package	PI	Principal Investigator
FCF	Fluids and Combustion Facility	PIV	Particle Image Velocimetry
FCU	FOMA Control Unit	PTCU	Payload Training Center Unit
FDSS	Fire Detection and Suppression System	REU	Remote Electronic Unit
FEA	Fluids Experiment Assembly	RFCA	Rack Flow Control Assembly
FIR	Fluids Integrated Rack	RHA	Rack Handling Adapter
FOMA	Fuel/ Oxidizer Management Assembly	RIP	Rack Interface Panel
FSAP	Fluid Science Avionics Package	RMSA	Rack Maintenance Switch Assembly
GIS	Gas Interface System	RPC	Remote Power Controller
GIU	Ground Integration Unit	RPCM	Remote Power Controller Module
GSE	Ground Support Equipment	RUP	Rack Utility Panel
HFR	High Frame Rate	S/ W	Software
HiBMs	High Bit depth Multi Spectral Camera Package	SAMS	Space Acceleration Measurement System
HR	High Resolution	SAR	Shared Accommodations Rack
H/W	Hardware	SSC	Station Support Computer
IAM	Image Acquisition Module	UIP	Utility Interface Panel
IOP	Input/ Output Processor	UV	Ultraviolet
IPP	Image Processing Package	VVS	Vacuum Vent System
IPSU	Image Processing and Storage Unit	WFCA	Water Flow Control Assembly
IR	Infrared	WTCS	Water Thermal Control System



# ***Fluids and Combustion Facility***

## ***Preliminary Design Review***



### **FCF Mission**

Support sustained, systematic microgravity fluid physics and combustion research on board the ISS in the following areas:

#### **Fluid Physics**

- Thermocapillarity
- Interfacial Phenomena
- Phase Change
- Complex Fluids
- Colloids
- Morphology
- Electrohydrodynamics
- Diffuso-Capillary
- Thermo-Diffuso-Capillary
- Multiphase Flow
- Granular Media

#### **Combustion**

- Laminar Flames
- Turbulent Combustion
- Reaction Kinematics
- Condensed Phase Organic Fuel Combustion
- Flame Spread and Fire Suppressants
- Droplet and Spray Combustion
- Soot and Polycyclic Aromatic Hydrocarbons



# ***Fluids and Combustion Facility***

## ***Preliminary Design Review***



### **FCF Overview**

- The FCF is a modular, multi-user facility being developed by GRC for permanent installation in the ISS US Lab Module to support microgravity fluid physics and combustion science.
- FCF establishes an initial operational capability for microgravity fluid physics and combustion science on-board ISS. The Project's scope is through the design, development, production, deployment and initial operation in ISS.
- The FCF Project transitions to a steady-state utilization program after on-orbit acceptance of the last FCF rack deployed.
- The FCF system consists of a Flight Segment and a Ground Segment. The Flight Segment consists of three powered racks and sufficient on-orbit stowage in ISS for FCF and payload hardware.
- FCF operates together with payload experiment equipment, ground-based operations facilities and the FCF ground segment to perform fluids and combustion science experiments within available FCF and ISS resources.
- The FCF Ground Segment is required for operation of the FCF and for successful integration and operation of experiment hardware and software into the FCF Flight Segment.





# Fluids and Combustion Facility

## Preliminary Design Review



### Incremental Deployment of the FCF System



Conduct Combustion Experiments only on orbit.



Conduct Combustion Experiments and Fluids Experiments on orbit.



Conduct Combustion, Fluids Experiments And complete flight segment requirements.




# Fluids and Combustion Facility

## Preliminary Design Review



### Three-Tier Solution for a Diverse Complement of Science Requirements

#### Fluids and Combustion Facility




#### Tier 1

Systems commonly needed by nearly all Fluids and Combustion Experiments

- Build once, launch once, use forever

#### Multi-Use Experiment Hardware




#### Tier 2

Multi-use experiment inserts customize FCF to a specific sub-discipline

- Build once, launch once per ~4 expts, reuse
- Other sub-discipline inserts reuse some equipment and concepts

#### PI Unique Equipment



#### Tier 3

PI unique equipment

- Build and launch for each experiment
- May reuse or add to capability of FCF or multi-use inserts




# Fluids and Combustion Facility

## Preliminary Design Review



### Three-Tier Solution for Combustion



#### PI-Specific Equipment

- Fuels/Oxidizers/Consumables
  - Gaseous Fuels/Oxidizers
  - Liquid/Solid Fuel Supplies
  - GC check gases
  - Adsorber cartridge filters
- PI-Unique H/W & S/W
  - Specific avionics/software
  - Specialized diagnostics or modules (e.g., cameras, lenses, filters, optics, control avionics and software)
- Specific Instrumentation

#### Multi-Use Insert (e.g., MDCA)

- Combustion Experiment Insert
  - Experiment orientation
  - Support structure
  - Flow ducts; test sections
- PI Avionics Package
  - Multi-use control electronics, motor drivers, etc.
  - Software modules
- Multi-Use Diagnostics / Lighting
  - e.g., Radiometers
  - PI-Unique illumination sources

#### FCF Combustion Integrated Rack

- Combustion Chamber
- Diagnostics
  - 8 standard locations (w/ SAR)
  - Image Processing & Storage
- Gaseous Fuel/Oxidizer Mgmt.
  - Gas Supply & Blending
  - Exhaust Vent / GC
- Electrical Power Control Unit
  - Power Conditioning & Control
- Command & Data Handling
- Active Rack Isolation
- Environmental Control
  - Water/Air Cooling
  - Gas Interface/Fire Detection




# Fluids and Combustion Facility

## Preliminary Design Review



### Three-Tier Solution for Fluid Physics



#### PI Specific Samples

- Samples with supporting hardware
- Specific Conditioning
- Specific Diagnostic

#### Investigation Module (e.g., LMM)

- Science Infrastructure (hardware/software) items that uniquely meet the needs of the PIs
- Unique Diagnostics
- Specialized Imaging
- Fluid Containment

#### Fluids Integrated Rack

- Power Supply
- Avionics/Control
- Common Illumination
- PI Integration Optics Bench
- Imaging and Frame Capture
- Fluid Diagnostics
- Environmental Control
- Data Processing
- Frangibles and Laser Light Containment






# Fluids and Combustion Facility

## Preliminary Design Review



### FCF Shared Accommodation Rack






# ***Fluids and Combustion Facility***

## ***Preliminary Design Review***




### **FCF Flight and Ground Segments**




#### **Flight Unit**

FCF racks located on-orbit in the ISS US Laboratory Module.



#### **Experiment Development Unit (EDU)**

High fidelity FCF model located at GRC. Available to payload developers for a variety of tests including interface verification, preliminary configuration selection, test sequence determination. Also used for sub-rack payload engineering model testing.



#### **Ground Integration Unit (GIU)**

FCF flight-equivalent racks located at GRC. The GIU is used for final flight interface verification between sub-rack payloads and the FCF, as well as acceptance testing of PI hardware. The GIU configuration may be maintained for on-orbit troubleshooting. The GIU is flight equivalent, except aluminum ISPRs are used and ARIS is not functional in the GIU.



#### **Payload Training Center Unit (PTCU)**

FCF training racks are located at JSC and used to train astronauts on FCF and sub-rack payload flight procedures. PTCU racks are supplemented with partial task trainers to prepare the ISS crew for experiment-specific operations.




# Fluids and Combustion Facility

## Preliminary Design Review



### FCF Ground and Flight Segments

#### ISS Fluids and Combustion Facility



#### Flight Segment


##### Space Station Based Equipment

- Combustion Integrated Rack
- Fluids Integrated Rack
- Shared Accommodations Rack
- Spares & Other Equipment

#### Ground Segment

##### Earth Based Operations

- Ground Integration Unit
- Experiment Development Unit
- PTC Training Unit
- Ground Handling & Testing Equipment
- FCF-Unique Telescience Equipment
- Operations Phase Processes





# ***Fluids and Combustion Facility Preliminary Design Review***

---



## **FCF System**

- Development Approach
- Design Drivers
- Design Philosophy
- System Design Features
- System Description







# Fluids and Combustion Facility

## Preliminary Design Review



### FCF Development Approach





# ***Fluids and Combustion Facility***

## ***Preliminary Design Review***

---



### **Key constraints that drive the FCF design**

- Rack Mass – limited by launch and on-orbit requirements
- Volume – limited to three contiguous racks plus one rack for stowage
- Power – limited by project requirements and ISS power resources
- Energy - limited by project requirements and ISS power resources
- Crew Time – limited amount of crew time and overall time to perform the FCF mission
- Number of Microgravity days
- Upmass - limited number of flights to ISS and limited upmass for PI and FCF hardware on each flight
- Data Communication – limited bandwidth for downlink of experiment data
- Stowage - limited stowage for on-orbit spares



# ***Fluids and Combustion Facility***

## ***Preliminary Design Review***



### **FCF Specification Key Requirements**

- Throughput - 5 combustion and 5 fluids typical basis experiments per year within the following resource constraints:
  - Power - 2000 W
  - Energy - 9000 kW-hrs
  - Thermal Cooling Water - 90 kg/hr
  - Data - 6 Mb/sec, 20 GB/day downlink
  - Crew time - 145 hrs/year
- Provide microgravity environment per experiment requirements
- Average Cost/PI - \$ 3.3M in 2005 US dollars
- 10-15 year on-orbit life
- Accommodation of first 6 real experiments
- Ability of three racks to operate as a three rack system

**Delivering the required science capability while working within the many constraints and requirements levied on the FCF is accomplished using rigorous systems engineering in a Concurrent Engineering environment.**




# Fluids and Combustion Facility

## Preliminary Design Review



### FCF Concurrent Engineering Approach





# ***Fluids and Combustion Facility Preliminary Design Review***

---



## **FCF Design Philosophy**

- 1) Design for highly automated operation
- 2) Maximize commonality across all three racks
- 3) Implement modular design
- 4) Employ design features for ease of operation, reconfiguration, maintenance and upgrade
- 5) Maximize use of COTS hardware and software wherever practical



# ***Fluids and Combustion Facility***

## ***Preliminary Design Review***

---



### **FCF Design Philosophy**

#### 1) Design for highly automated operation

FCF hardware and software is designed for automated operation from the ground once the crew has completed the experiment set up.

- Provides greater flexibility for science operations by the PI
  - revision of experiment timeline to react to unexpected science phenomena
  - rerun test points that require a second look
- Dramatic reduction in crew time required to complete experiment test matrix



# ***Fluids and Combustion Facility Preliminary Design Review***



## **FCF Design Philosophy**

### 2) Maximize commonality across all three racks

Maximizing commonality across all three racks and between subsystems reduces program risk while optimizing for the fluids and combustion missions.

- Reduces common hardware, software, test methods, inspections, assembly and test procedures, materials, etc. improves repeatability and product quality
- Common operational procedures reduce risk of crew confusion
- Reduces total number of analyses, tests, and designs
- Reduces total life cycle costs
- Allows common sparing on orbit




# Fluids and Combustion Facility

## Preliminary Design Review



### FCF Common Hardware


#### Structural




Slides




Rack Door



Left and Right  
Pin Assembly




Center  
Post




Removable Latch Mount


#### Command and Data




IPSU



IOP




IAM




DCM


#### GFE



EPCU




ARIS and ISPR



SAMS FF

#### ECS



ATCU with ATCU Rack  
Attachment Hardware

GIS

WTCS

FDSS





# ***Fluids and Combustion Facility***

## ***Preliminary Design Review***

---



### **FCF Design Philosophy**

#### 2) Implement modular design

Modular design of hardware and software:

- Allows interchange of Orbital Replacement Units (ORUs) between racks, increasing system flexibility and reducing stowage needs
- Reduces crew time
- Aids in sparing
- Provides for fault isolation
- Accommodates new ORU upgrades or new PI hardware and software
- Modular software approach uses latest development tools/techniques, object oriented analysis, and design philosophy to maximize commonality, reusability and maintainability.



# ***Fluids and Combustion Facility***

## ***Preliminary Design Review***

---



### **FCF Design Philosophy**

- 3) Employ design features for ease of operation, reconfiguration, maintenance and upgrade. Examples include:
- 1/4 turn captive fasteners
  - quick disconnects for fluid lines
  - C-mount adapters for diagnostic packages
  - one handed package change out using quick latch device
  - waffle plate bench with quick thumb latch hold down in FIR and SAR




# ***Fluids and Combustion Facility***

## ***Preliminary Design Review***



**Modular design of packages along with quick latch device provides fast, easy reconfiguration of the bench.**





# ***Fluids and Combustion Facility***

## ***Preliminary Design Review***

---



### **FCF Design Philosophy**

4) Maximize use of COTS hardware and software wherever practical to reduce life cycle costs and provide readily available upgrade paths.

Examples include:

- processor boards and disk drives
- diagnostic elements (cameras, lenses, etc.)
- VxWorks operating system
- fluid supply bottles
- gas chromatograph
- others



# ***Fluids and Combustion Facility***


## ***Preliminary Design Review***



### **FCF Design Feature: Deployable Optics Bench**

Deployable optics bench:

- Allows crew easy access to both front and back of bench
- Minimizes time required to reconfigure rack between experiments






# ***Fluids and Combustion Facility***

## ***Preliminary Design Review***




### **FCF Design Feature: Commonality between Optics Bench designs**

Cooling, Power, Data and Mechanical Interfaces allow for interchangeability of packages between racks.



CIR Optics Bench



FIR & SAR Optics Bench



# Fluids and Combustion Facility


## Preliminary Design Review




### FCF Design Feature: Ducted Bench Cooling

Ducted bench cooling:

- Provides ease for PI design from a thermal standpoint
- Better control of ATCU exit air
- Reduces number of fans and other components that are potential noise and microgravity disturbance generators



CIR uses positive pressure



FIR and SAR use negative pressure





# ***Fluids and Combustion Facility***


## ***Preliminary Design Review***



### **FCF Design Feature: No Tools Required for Reconfiguration**

Standard connections for components and packages results in reconfiguration between experiments:

- Simplifies crew procedures
- Enhances ability to use between racks
- Minimizes stowage
- Reduces crew time required for reconfiguration between experiments



**IAM Four Flange Mount Interface**





# Fluids and Combustion Facility


## Preliminary Design Review



### FCF Design Feature: Rack to Rack Communication

Communication between racks:

- Allows better use of resources
  - Offloading of mass
  - Space for additional diagnostics
- Allows for centralized control of all three racks






# Fluids and Combustion Facility

## Preliminary Design Review



### FCF Design Feature: Rack Electrical Grounding



- ISS Secondary power isolated from chassis
- Separately derived electrical power source connected to structure at no more than one point
- Package-to-package communication designed to eliminate ground loops utilizing fiber optics or isolated differential signals
- Wiring and cabling shall meet applicable requirements specified in SSP 30242, ISS Cable/Wire Design and Control Requirements for Electromagnetic Compatibility
- System grounding shall meet applicable requirements specified in SSP 30240, ISS Grounding Requirements




# Fluids and Combustion Facility

## Preliminary Design Review



### FCF Design Feature: Rack Bonding



- Entire optics bench mounting surface shall be electroless nickel.
- Diagnostic package mounting surface to optics bench shall be electroless nickel.
- ARINC optics bench and diagnostic package connectors surface material shall be electroless nickel.
- Bonding strap shall be used to establish an electrical bond between the optics bench and rack per SSP 30245.
- Bonding straps shall be attached to the nickel-plated mounting locations provided on the rack corner post.
- On-orbit bond reverification of interchangeable optics bench packages will be accomplished via a visual inspection of the mating surfaces. Bonding is accomplished through mounting surfaces of the diagnostics package and the optics bench. A secondary path is through the connector shell.
- Packages separate from the optics bench such as the EPCU and IOP will utilize bonding straps attached to the rack corner post.




# *Fluids and Combustion Facility*

## *Preliminary Design Review*



### **FCF Design Feature: Fiber Optic Serial Data Link (SDL)**

- SDL converts parallel data lines from diagnostic packages onto a very-high speed optical transmission medium. At the receiving end the data is restored to a parallel format for processing.
- Using the SDL greatly reduces the number of conductors in each rack and allows rack to rack communications.
- Fiber optic line is very immune to any noise interference and does not radiate any electrical noise.
- PI provided digital cameras may require customization of SDL at the camera interface.





# ***Fluids and Combustion Facility***

## ***Preliminary Design Review***

---



### **FCF System**

#### **Each rack has stand alone capabilities**

- CIR will be able to conduct combustion experiments.
- FIR will be able to conduct fluids experiments.
- SAR will be able to conduct fluids experiments, combustion experiments (if required containment vessel is provided by PI), and middeck experiments.

#### **The integrated system has enhanced capabilities**

- Augment FIR capabilities with SAR in order to conduct range of fluids experiments that cannot be fully accomplished by the FIR alone.
- Augment CIR capabilities with SAR in order to conduct range of combustion experiments that cannot be fully accomplished by the CIR alone.
- Improve FIR and CIR efficiency by SAR processing data and down linking previous run while experiments continue in FIR or CIR. This will allow for more science to be performed.



# ***Fluids and Combustion Facility***

## ***Preliminary Design Review***

---



### **FCF System**

#### **Benefits of the Integrated FCF system**

#### **SAR completes the FCF system**

- Provides flexibility to off-load hardware from optics plates in CIR and FIR to create additional space for PI-specific equipment, enhancing science data return.
  - CIR - opportunity to view into chamber from all 8 window locations
  - FIR - enables simultaneous operation of high-resolution, high frame rate diagnostics
- Capability to provide additional data storage in the SAR allows for greater number of test runs that would otherwise be precluded by ISS downlink limitations and/or a single rack's data storage capability.
- Allows for experiments that would otherwise exceed a single rack's limit for mass, volume or cooling capability.






# Fluids and Combustion Facility

## Preliminary Design Review




The SAR Permits FCF to Operate at Full Performance



### Advanced FCF Data Handling - *Completes FCF Imaging/Data Capabilities*

- **Input/Output Processor** for facility command, control, health and status monitoring
- **Image Processing Storage Units (2)** for combustion and fluids image acquisition, processing, storage
- **Additional Data Processing Capability** for processing data from multiple experiments
- **Mass Storage** of FCF Science Data with removable media for return to Earth
- **Video Monitor** for crew interface
- **SAR Science Avionics Package** for experiment data handling.
- **Imaging Interface** for other US Lab science rack usage.



### Additional FCF Science Accommodations - *Science Support and Throughput*

- **Science Optics Bench** for hardware mounting, servicing and access, as well as experiment accommodation
- **Science Sample Conditioning** (e.g., powered fluid sample preparation and storage)
- **Mid-deck Locker Interface** with avionics simulators and interfaces

### Potential for FCF Upgrades — *Enhanced Capabilities Possible after SAR Deployment*

- **Calibration Apparatus** for FCF equipment requiring regular on-orbit calibration.
- **On-Orbit Logistical Equipment** to maximize system readiness & throughput
- **FCF Diagnostic Compatibility** with CIR and FIR diagnostics to support SAR science payloads
- **Multiple PI-Specific Hardware** available on-orbit at all times to maximize flexibility and thru-put
- **Gas/Fluid Supply** to support CIR/FIR science
- **Chemical Analysis Capability** for on-orbit fluid sampling and analysis (e.g., mass spectrometer)

The FCF Shared Accommodations Rack (SAR) provides centralized command & data handling capability and multipurpose science volume needed to maximize FCF system effectiveness and meet fluids and combustion science requirements.




# Fluids and Combustion Facility

## Preliminary Design Review



### FCF System








# Fluids and Combustion Facility

## Preliminary Design Review



### FCF System






# Fluids and Combustion Facility

## Preliminary Design Review



### FCF System: Interfaces

Functional ISS  
To FCF to  
Payload  
Diagram  
(flight segment)



\* - including analog & digital video




# Fluids and Combustion Facility

## Preliminary Design Review



### FCF System: Interfaces

#### ISS to Rack Interfaces through UIP






# Fluids and Combustion Facility

## Preliminary Design Review



### FCF System: Interfaces

#### FCF Electrical System: Rack to Rack Interface





# Fluids and Combustion Facility

## Preliminary Design Review




### FCF System: Interfaces

FCF Electrical System: Rack to Rack Interface

#### Rack to Rack Umbilical Concept

- ✓ Fiber optic cable bundle routed from SAR to each science rack along z-panel.
- ✓ Connectors added to ARIS pass through panel at each rack.
- ✓ FCF diagnostic capabilities are accessible to non-FCF payloads via connectors on SAR, CIR, FIR pass-through panels.
- ✓ Allows high-speed data transfer between racks independent of ISS C & DH system.





# ***Fluids and Combustion Facility Preliminary Design Review***

---



## **Analysis Approach and Processes**





# ***Fluids and Combustion Facility Preliminary Design Review***

---



## **Analysis and Test Approach**

- Major system approaches have been described for each of the major flight environments in plans that integrate the analysis and testing for the following critical areas:
  - Microgravity Control
  - Acoustic Noise Control
  - Thermal Analysis and Test
  - Structural Analysis and Vibration Testing
  - EMI Compatibility




# Fluids and Combustion Facility

## Preliminary Design Review



### Microgravity Control Approach






# Fluids and Combustion Facility

## Preliminary Design Review



### Microgravity Verification – Current Status






# ***Fluids and Combustion Facility Preliminary Design Review***



# Initial Process

# Acoustic Noise Control Approach






# Fluids and Combustion Facility

## Preliminary Design Review



### Acoustic Noise Verification - Current Status







# Fluids and Combustion Facility

## Preliminary Design Review



### Thermal Analysis Overview








# Fluids and Combustion Facility

## Preliminary Design Review




### Thermal test approach





# Fluids and Combustion Facility

## Preliminary Design Review



## FCF STRUCTURAL ANALYSIS PROCESS



# ***Fluids and Combustion Facility***

## ***Preliminary Design Review***

---



### **STRUCTURAL ANALYSIS APPROACH**

- Flight safety requires analysis of all safety critical structures
- Safety critical structures defined as structure in payload primary load path including pressurized systems
- Requirement met using a combination of classical (hand) and finite element analysis
  - Classical analysis used on fasteners (Chambers approach) and certain parts like FOMA manifolds
  - Finite element models created for each FCF subsystem (Optics Bench Assembly, ATCU, etc.)
- Fracture critical parts assessed using FLAGRO
- Fatigue life assessed using fracture analysis (FLAGRO) – a conservative approach



# ***Fluids and Combustion Facility Preliminary Design Review***



## **STRUCTURAL ANALYSIS APPROACH (cont.)**

- Finite element models integrated into ISPR model for rack level analysis
  - Launch and landing loads stress analysis
  - Launch and landing coupled loads analysis
  - ARIS performance analysis
  - Microgravity disturbance analysis
- Few separate (less detailed) dynamic FE models; stress model fidelity required for higher frequency analysis
- For PDR a summary of results are documented in Safety Critical Structures Data Package
- For CDR a detailed stress report will be written for each major structural element




# Fluids and Combustion Facility

## Preliminary Design Review



### FCF Approach to Random Vibration Testing Flight and Qualification Hardware






# Fluids and Combustion Facility

## Preliminary Design Review



Safety Critical/Fracture Critical Part Screening Process  
for Non-rotating, Metallic Components








# Fluids and Combustion Facility

## Preliminary Design Review



### Safety Critical/Fracture Critical Selection Logic





# ***Fluids and Combustion Facility***

## ***Preliminary Design Review***

---



### **EMI Control Approach**

- FCF Electromagnetic Interference/Compatibility Plan (FCF-PLN-0027) has been developed to provide requirements for EMI/EMC testing at both the subsystem and rack level.
- EMC compliance will be accomplished in two phases.
  1. Subsystem Level – EPCU as GFE, IOP, FCU IPP, DCM, PI, etc.

Subsystem level testing may also include preliminary tests at the breadboard level for hardware believed to be at risk.
  2. Individual Rack Level – CIR, FIR, and SAR

Rack level requirements are well defined and understood. They have been performed on other ISS hardware by the GRC EMI Laboratory. Test levels and limits are defined in SSP 57000 and SSP 30237. Test methods are defined in SSP 30238.
- This subsystem design approach will be verified via successful rack level tests.
- Subsequent testing at the subsystem or box level will be used as justification for on-orbit replacement of existing hardware or the addition of new PI hardware.
- In cases where the hardware does not pass subsystem level test requirements, an analysis of the out of specification condition will be performed. This analysis will assess the potential impact on other subsystems as well as expected performance at the rack level.



# Fluids and Combustion Facility

## Preliminary Design Review



### Status of FCF Analyses

- Initial Analysis
  - Functional (SRED compliance)
  - 57000 requirements
  - Safety Requirements
  - “A” Specification
- Preliminary Analysis
  - Functional (SRED compliance)
  - Structural (at the rack level)
  - Thermal (at the rack level)
  - Acoustic Noise
  - Microgravity including ARIS effects
  - Fit check of all components using Pro-E
  - 57000 requirements

Complete
Complete
Complete
Complete
Complete
In Work
In Work
In Work
In Work
Complete
Complete



# ***Fluids and Combustion Facility Preliminary Design Review***

---



## **Development/Assembly and Integrated Test Approach and Processes**




# Fluids and Combustion Facility

## Preliminary Design Review




### Development, Assembly and Integration Approach





# Fluids and Combustion Facility

## Preliminary Design Review









# Fluids and Combustion Facility

## Preliminary Design Review



### EM Development & Test

#### Selective Functional Testing



- Each package undergoes qualification and acceptance testing to verify that its design meets required performance parameters.
- Qualification and acceptance level testing of **all** packages accomplished using EM hardware configured in one of three ways, as appropriate:
  - as stand alone package
  - as part of a larger subsystem
  - as part of a partially integrated rack
- Each package and fully integrated rack also has a full range of EM acceptance testing.
  - EMI, acoustics,  $\mu\text{g}$ , power draw, pressure, thermal performance



# Fluids and Combustion Facility

## Preliminary Design Review



### EM Development & Test - Qualification Testing

- Qualification Thermal Cycling testing is performed at the stand alone package level.
- Qualification vibration testing is performed with the hardware in an appropriate test configuration to ensure that test conditions reflect the real vibration environment the packages will see during launch. For the purpose of example, CIR vibration test configurations are shown.

#### Stand Alone Package

Unstowed at launch

- EPCU
- IOP

Stowed at launch (a)

- IPP
- IPSU
- GC
- FCU
- Diagnostic Pkgs

#### Integrated Subsystem

- Optics Bench, FOMA Pkgs & Chamber (CIR)\*

*\* All packages mounted on Optics Bench in the launch configuration*

#### Partially Integrated Rack

- ECS Components, FDSS, GIS and Rack Doors

#### Fully Integrated Rack


- No vibration qualification testing at the fully integrated rack level

(a) Packages that are stowed during launch will be qual tested in their launch configuration; i.e., wrapped in foam and placed in a stowage container.



# Fluids and Combustion Facility


## Preliminary Design Review





# Fluids and Combustion Facility

## Preliminary Design Review



Launch Configuration




Stowed Hardware



# Fluids and Combustion Facility

## Preliminary Design Review



Launch Configuration



Stowed Hardware




# Fluids and Combustion Facility

## Preliminary Design Review



Use of PRCU and GIU to interface with flight racks for integrated testing





# ***Fluids and Combustion Facility Preliminary Design Review***

---




## **Fluids Combustion Facility Flight and Ground Safety**





# Fluids and Combustion Facility

## Preliminary Design Review





# ***Fluids and Combustion Facility Preliminary Design Review***



## **Fluid Combustion Facility (FCF) Flight Safety Review/Package Approach Combustion Integrated Rack (CIR)**

- Safety Reviews Planned
  - Flight - Phase 0/1, II Reviews will be held independently
    - Phase III Flight Review will integrate the MDCA Facility Package
- Flight - Phase 0/1, II Safety Data Packages(SDP) will be prepared independently
  - Phase III Flight SDP will incorporate MDCA Phase III SDP

## **Fluids Integrated Rack (FIR)**

- Safety Reviews Planned
  - Flight - Phase 0/1, II Reviews will be held independently
    - Phase III Flight Review will integrate the LMM Facility Package
- Flight - Phase 0/1, II Safety Data Packages(SDP) will be prepared independently
  - Phase III Flight SDP will incorporate LMM Phase III SDP

## **Shared Accommodation Rack (SAR)**

- Safety Reviews Planned
  - Flight - Phase 0/1, II Reviews will be held independently
    - Phase III Flight Review will integrate the SAR Facility Package
- Flight - Phase 0/1, II Safety Data Packages(SDP) will be prepared independently
  - Phase III Flight SDP will incorporate SAR Phase III SDP



# ***Fluids and Combustion Facility***

## ***Preliminary Design Review***

---



### **FCF Payload Flight Safety Review/Package Approach**

#### **Multi-User Droplet Combustion Apparatus (MDCA)**

- Safety Reviews Planned
  - Flight - Phase 0/1, II Reviews will be held independently
    - Phase III Flight Review will be integrated with CIR Facility
- Flight - Phase 0/1, II Safety Data Packages(SDP) will be prepared independently
  - Phase III Flight SDP will be prepared and integrated with CIR Phase III SDP

#### **Light Microscopy Module (LMM)**

- Safety Reviews Planned
  - Flight - Phase 0/1, II Reviews will be held independently
    - Phase III Flight Review will be integrated with FIR Facility
- Flight - Phase 0/1, II Safety Data Packages(SDP) will be prepared independently
  - Phase III Flight SDP will be prepared and integrated with FIR Phase III SDP

#### **Payload for Shared Accommodations Rack (SAR)**

- Safety Reviews Planned
  - Flight - Phase 0/1, II Reviews will be held independently
    - Phase III Flight Review will be integrated with CIR Facility
- Flight - Phase 0/1, II Safety Data Packages(SDP) will be prepared independently
  - Phase III Flight SDP will be prepared and integrated with SAR Phase III SDP



# ***Fluids and Combustion Facility Preliminary Design Review***



## **Fluids Combustion Facility/Payload Ground Safety Approach**

- Multi-User Droplet Combustion Apparatus (MDCA) Safety Reviews Planned
  - Ground - Phase II & III
    - Reviews will be held in conjunction with the CIR Facility Ground Review
  - Ground Payload Safety Review Packages will be prepared Independently
- Light Microscopy Module (LMM) Safety Reviews Planned
  - Ground - Phase II & III
    - Reviews will be held in conjunction with the FIR Facility Ground Review
  - Ground Payload Safety Review Packages will be prepared Independently
- Payload TBD for Shared Accommodations Rack (SAR) Safety Reviews Planned
  - Ground - Phase II & III
    - Reviews will be held in conjunction with the SAR Facility Ground Review
  - Ground Payload Safety Review Packages will be prepared Independently



# ***Fluids and Combustion Facility***

## ***Preliminary Design Review***

---



### Flight Safety Status

- Flight Facility Safety
  - Combustion Integrated Rack (CIR) – Flight Safety
    - CIR Phase 0/1/Delta Phase 0/1 Flight Safety Review Completed Successfully
    - CIR Phase II Flight Safety Review will be scheduled for TBD
  - Fluids Integrated Rack (FIR) – Flight Safety
    - FIR Phase 0 Flight Safety Review completed successfully.
    - FIR Phase I Flight Safety Review completed successfully.
- Flight FCF Payload Safety
  - Multi-User Droplet Combustion Apparatus (MDCA) – Flight Safety
    - MDCA Phase 0/1 Flight Safety Review Completed
    - MDCA Phase II Safety Review will be scheduled for TBD
  - Light Microscopy Module (LMM) – Flight Safety
    - LMM Phase 0/1 Flight Safety Review scheduled for March, 2001
    - LMM Phase II Flight Safety Review scheduled for March, 2001



# ***Fluids and Combustion Facility Preliminary Design Review***



## **Ground Safety Status**

### **Ground Facility Safety**

- Combustion Integrated Rack (CIR) – Ground Safety
  - CIR Phase II Ground Review Scheduled for TBD
- Fluids Integrated Rack (FIR) – Ground Safety
  - FIR Phase II Ground Review Scheduled for TBD

### **Ground FCF Payload Safety**

- Multi-User Droplet Combustion Apparatus (MDCA) – Ground Safety
  - MDCA Phase Draft Ground Safety Package Completed Review to be scheduled.
  - MDCA Phase II Ground Review Scheduled for TBD
- Light Microscopy Module (LMM) – Ground Safety
  - LMM Phase Draft Ground Safety Package to be completed.
  - LMM Phase II Ground Review Scheduled for TBD



# ***Fluids and Combustion Facility Preliminary Design Review***

---



## Safety Key Requirements Compliance

- No non-compliances
- FIR Phase I Review - January 30 & 31, 2001
  - System Laser Safety approach
  - Fluids containment supplied by PI
- CIR Phase 0/1 complete, Phase II in July 01
  - FOMA system which provides gases has been through safety process and White Sands O2 analysis
  - CIR provides Chamber and equivalent 2 levels of containment
  - System Laser Safety approach





# ***Fluids and Combustion Facility Preliminary Design Review***

---



## **Differences Between EDU, GIU and Flight Segment**



# ***Fluids and Combustion Facility Preliminary Design Review***



Differences between EDU, GIU and Flight Segment

- Breadboard

- Hardware

- New hardware and EM hardware

- Hardware functionally similar to flight hardware

- Supported with necessary GSE (including structural support)

- Function

- Will functionally similar to flight unit

- Used by PI's to test feasibility of hardware/concepts System Specification



# ***Fluids and Combustion Facility Preliminary Design Review***



## Differences between EDU, GIU and Flight Segment

- EDU

- Hardware

- Qualification units, Functional equivalent and EM hardware

- Complete set of hardware as in flight unit (sans ARIS and SAMS)

- Supported with necessary GSE (including structural support)

- Function

- Will functionally equivalent to flight unit

- Used by PI's to develop flight hardware



# ***Fluids and Combustion Facility***

## ***Preliminary Design Review***



### Differences between EDU, GIU and Flight Segment

- GIU

- Hardware

Flight Pedigree hardware. (Drawings will be flight hardware drawings).  
Hardware will be functionally tested as flight hardware was but not environmentally tested.

Quality will be same as flight unit.

Complete set of hardware as in flight unit (sans ARIS and SAMS)

Supported with necessary GSE (including structural support). PRCU will supply water, power and vacuum

- Function


Will functionally equivalent to flight unit.

Used by PI's for integration verification, Flight S/W update validation, trouble shooting, and following on orbit progress. Procedure development, verification of system level requirements between racks (see chart)



# Fluids and Combustion Facility

## Preliminary Design Review





# ***Fluids and Combustion Facility Preliminary Design Review***

---



## **FCF Level Risks**



# Fluids and Combustion Facility

## Preliminary Design Review



### Risk Management

- All risks for the FCF are identified and assessed by the FDC team. These risks are then tracked in the MRDOC automated Risk Status and Tracking System. Risks are reviewed weekly by the PDT Leads, assessing the progress of risk mitigation activities.

Microsoft Access - [Active Risk Database]

FDC FEDERAL DATA CORPORATION Major Programs Sector

**Risk Form**

Hernandez Engineering Inc.

Project Name: Table4

ID: 108 Date Opened: 10/4/00

Risk Name: [Text Box]

Risk Identifier: [Text Box]

Assignee: [Text Box]

Risk Statement

Condition: [Text Box]

Consequences: [Text Box]

Updated: [Text Box]

Status: [Dropdown]

L.C. Stage: [Dropdown] Add

Classification: [Dropdown] Add

Action Type: [Dropdown]

Enter risk manually?: ☐

Risk Level: 1 - 5, 1 = Highest

Impact: [Dropdown]

Time Frame: [Dropdown]

Probability: [Dropdown]

Prioritize Risks

Closure Criteria: [Icon] Action: [Icon] Tracking: [Icon] Context: [Icon] Progress: [Icon] Closure Details: [Icon]

Description of the risk and its effects

NUM

Microsoft PowerPoint - [Pr...] Active Risk Database

5:27 PM





# ***Fluids and Combustion Facility Preliminary Design Review***



## **Top Risks for FCF System**

The following issues represent significant risk to the FCF's ability to meet mission objectives. An FCF mitigation plan has been developed for each risk area:

- Microgravity disturbances
  - addressed by Microgravity Control Plan
- Acoustic noise
  - addressed by Acoustic Noise Control Plan
- Radiation effects on electronics
  - mitigation and control plan in work

Risk items associated with a particular rack are addressed at the rack level.

In all, 45 active risks are currently being tracked by the FDC Team.



# Fluids and Combustion Facility

## Preliminary Design Review



- Developed a formal Risk Management Plan FCF-PLAN-0046 which establishes structured risk management processes.
- Established a risk database to track and control project risks.

<b>Identify</b>	Capturing statements of risk Capturing the context of a risk
<b>Analyze</b>	Evaluating attributes of risk Classifying risks Prioritizing risks
<b>Plan</b>	Assigning responsibility Determining approach Defining scope and actions
<b>Track</b>	Related data is acquired, compiled, analyzed, and reported <i>Note: Risks can be tracked individually or in sets</i>
<b>Control</b>	Decisions made based on the data presented in the tracking reports <i>Note: Risks can be controlled individually or in sets</i>

### Continuous Risk Management



**Communication and documentation shall be performed through out the continuous risk management process**



# Fluids and Combustion Facility

## Preliminary Design Review



### Impact Definitions

Loss of entire system Severe Down Time > 30 days	<b><u>Catastrophic</u></b>
Major system damage Significant down time 1 – 30 days	<b><u>Critical</u></b>
Minor system damage Minor down time < 1 day	<b><u>Marginal</u></b>
No system damage Insignificant down time < 1 hour	<b><u>Negligible</u></b>

### Risk Priority Matrix

From MRD-DOC-0006 (Risk Management Plan)

Severity Level	Likelihood of Occurrence			
	Probable	Occasional	Remote	Improbable
Catastrophic	1	1	2	3
Critical	1	2	3	4
Marginal	2	3	4	5
Negligible	3	4	5	5

**Key**

1	HIGH
2	MEDIUM HIGH
3	MEDIUM
4	MEDIUM LOW
5	LOW




# Fluids and Combustion Facility

## Preliminary Design Review



### Flowchart For Radiation Effects Analysis



\*Standard FMEA will be performed on all packages. These FMEA's look specifically for SEE



# ***Fluids and Combustion Facility Preliminary Design Review***

---



## **Radiation Effects Analysis - Subsystem Priority Listing**

- Primary
  - Common IOP (Tom Young/Mike Casciani)
  - Common ECS (Tracy Neff/John Siamidis)
  - FIR FSAP (Tom Young/Bart Gruber)
  - CIR FCU (Myron Wilkinson/Mary Domo)
- Secondary
  - FOMA & Combustion Chamber (TBD)
  - CIR IPP, Common IPSU (TBD)
  - Science Diagnostics/DCM (TBD)



# ***Fluids and Combustion Facility Preliminary Design Review***

---



## **Radiation Effects Mitigation Plan**

- Environmental profile has been established.
- Gather data (parts lists) for initial analysis of radiation effects.
- Determine probability of occurrence.
- Test/Analyze based upon above factors.
- Mitigate through acceptance or redesign.



# ***Fluids and Combustion Facility Preliminary Design Review***

---



## **FCF System Compliance**





# ***Fluids and Combustion Facility Preliminary Design Review***



## Common Hardware Non-compliances

- Captive Fasteners on Several Components
  - Components include:
    - Optics Bench Attachment Hardware
    - Door Launch Restraints
    - Rack Attachment Hardware for ATCU, EPCU, and IOP
  - Mitigation strategy includes:
    - Identifying and assessing candidate alternatives
    - If necessary, rationale for potential waiver is one-time of operation of these fasteners for rack setup or ORU failure
    - Resolution should occur by April 2001



# ***Fluids and Combustion Facility Preliminary Design Review***



## **Common Hardware Non-compliances**

- Diagnostics Package Pin ID, Connector Arrangement, and Ease of Disconnect
  - All affected by use of Airborn Connectors on several packages
  - Rationale for seeking waiver is infrequency of operation of these connectors and inability to obtain packing density with larger connectors that meet current 57000 requirements
  - Waiver to be submitted by April 2001



# ***Fluids and Combustion Facility Preliminary Design Review***



## **Common Hardware Non-compliances**

- **IPP Accessibility**

- Installation and removal of the IPP involves mating and demating an electrical connector – the UML ARINC connector
- When installed in the CIR, the IPP will be bolted to the optics bench
- Because of their close proximity, it may not be possible to install or remove the IPP without removing adjacent diagnostic packages
- If crew consensus is obtained at the crew review, will seek a waiver of this requirement by June 2001



# ***Fluids and Combustion Facility Preliminary Design Review***

---



## **Summary of Issues**



# ***Fluids and Combustion Facility Preliminary Design Review***

---



## **FCF System Issues**

The following items are significant issues relative to the FCF's ability to meet requirements:

### ISS Resource Issues

- FCF system spec power constraint of 2000 watts
- Data downlink available for experiment data

### FCF Design Issues

- Diagnostics and other bench mounted components cannot meet Kick Load requirement



# ***Fluids and Combustion Facility Preliminary Design Review***

---



**Notes**